

Date: 01/01/2016

التاريخ: 2016/01/01

Service Agreement includes PPM & Emergency Calls

اتفاقية خدمة شاملة الصيانة الدورية والزيارات الطارئة

First Party: Magrabi center RD

الطرف الأول : مركز مغربي طريق المدينة

Second Party: Al-Amin medical Instruments Co. Ltd.
Country: Kingdom of Saudi Arabia - City: Jeddah
Telephone: 966-12-6601149
Fax: 966-12-6601146

الطرف الثاني : شركة الأمين للتجهيزات الطبية والعلمية
الدولة المملكة العربية السعودية - المدينة: جدة
هاتف : 966-12-6601149
فاكس : 966-12-6601146

With the help of Allah, the two parties have agreed on a maintenance contract for medical equipment on the following terms:

انه وبسوء الله تعالى اتفق الطرفان على عقد الصيانة للمعدات الطبية الموضحة ادناه و ذلك طبقا للشروط الآتية:

1. Description of work and equipment:

- The Second Party undertakes to carry out all basic maintenance works (Provided that the equipment is in 100% good operational condition before Implementation of contract). This Contract covers unlimited emergency calls. (Spare parts are excluded)

Accessories and consumables are not included in any type of our service Agreements.

- Location:
Jeddah, Makkah – K.S.A
- Equipment covered by contract:
Listed on Supplement 'A'

- وصف العمل والمعدات:-
 - يتعهد الطرف الثاني بالقيام بكافة اعمال الصيانة الأساسية (بشرط أن تكون الأجهزة بحالة تشغيلية 100% قبل تنفيذ هذه الاتفاقية) التي تشمل أجور العمالة للصيانة الدورية وعدد غير محدود من الزيارات الطارئة (ولا تتضمن قطع الغيار)

علما بان الاكسسوارات والمستهلكات غير مشمولة بأي نوع من عقود صيانتنا.

- الموقع:
جدة - المملكة العربية السعودية
- المعدات موضوع العقد:
متضمنة في ملحق أ

2. Contract term: Qty "1" Gregorian year
Starts on: 01/01/2016
Ends on: 31/12/2016

2- مدة العقد:- عدد "1" عام ميلادي
تبدأ في : 2016/01/01م
تنتهي في : 2016/12/31م

3- # of Visits: "2" visits / Year

3- عدد الزيارات : "2" زيارة سنويا

4. Contract amount: 150 000 Saudi riyals

4- قيمة العقد: 150 000 ريال سعودي

5. Payment method:
Half the value of the contract upon sign the contract.
Half the value of the contract after six months in 01/06/2016

5- نظام الدفع:-
نصف قيمة العقد عند توقيع العقد
نصف قيمة العقد بعد ستة أشهر في 2016/06/01.

6. Preventive maintenance
The Second Party shall undertake all preventive, remedial maintenance and calibration for equipment mentioned above as per the PPM schedule in supplement "A" in accordance with manufacturer's recommendations as per attached check list In supplement "B" using qualified Engineers and Technicians of the best quality and experience.

6- الصيانة الوقائية :-
يقوم الطرف الثاني باجراء جميع اعمال الصيانة الوقائية والتصحيحية والمعايرة للجهاز المذكور أعلاه وذلك حسب المواعيد المذكورة بالملحق "A"، وطبقا للتوصيات الخاصة بالشركات المصنعة كما هو موضح بالملحق "B"، وبأعلى أداء وبواسطة مهندسين وفنيين على مستوى عالي من الكفاءة والخبرة.

7. Second Party shall guarantee level of rendered maintenance as follows:

Company engineer to respond to maintenance request within 48 work hours, latest effective receipt of request by Fax or e-mail from the customer's authorized officer.
Maintenance that does not require spare parts, or in case spare parts are available within the Kingdom, company shall guarantee completion of the required repair within 96 work hours or less from receipt of maintenance request from the authorized officer of

7- يضمن الطرف الثاني مستوى الصيانة المقدم كما يلي :-
سيقوم مهندس الشركة بالاستجابة لأي طلب صيانة خلال مدة أقصاها 48 ساعة عمل من استلام طلب صيانة سواء بالفاكس أو بالبريد من الشخص المخول بطلب الصيانة من المستشفى.
بالنسبة للصيانة التي لا تحتاج إلى قطع غيار أو في حالة توفر قطع الغيار المطلوبة في المملكة فإن الشركة تضمن اكمال الإصلاح المطلوب خلال 96 ساعة عمل أو اقل من استلام طلب صيانة سواء بالفاكس أو بالبريد من الشخص المخول بطلب الصيانة من الطرف الأول

the customer.

- For repairs that require spare parts not available in Kingdom, the company guarantees completion of repair within 8 workweeks from receipt of maintenance request by Fax or e-mail from the authorized officer of the customer.

Company undertakes to provide highly qualified trained engineers for equipment covered by this contract and sold by the Company to First Party.

8. Fault reporting system

First Party shall notify the Second Party of faults by Fax, E-mail or on phone, to the maintenance department of the company. Reports shall indicate the type and date of fault. Faults shall be reported during company work hours from 08:30AM to 01:30PM and from 02:30PM to 07:00 PM Exception of this will be Fridays, Saturdays and official holidays.

- Call Center No. : 920026426 .

9. Spare parts

The Second Party shall supply the spare parts required by the first party for repair of equipment included in this contract, based on a report to be submitted by the engineer of the Second Party and a quotation to be accepted by the First Party

10. Maintenance service

The Second Party shall provide the maintenance services of the equipment mentioned in this contract. Such services shall be provided during regular work hours during weekdays from Sunday to Thursday : from 08:30AM to 01:30PM and from 02:30PM to 07:00 PM Maintenance services provided by the company shall cover the following:

- Undertaking scheduled service inspections as specified by the manufacturer, provided that such inspections are performed as previously contained in agreement, including lubrication, cleaning, equipment function inspection, current leak inspection, setting and calibration of equipment. Inspection timing shall be agreed in advance in such a way as to provide company technicians' full access to the equipment listed in contract to undertake such inspections. In case company officials are denied such access to equipment, inspection work shall be suspended.

Exceptions:

This agreement shall not cover maintenance and repair services resulting from causes beyond the control of the company, such as the acts of God, acts of civil and military authorities, fires, civil commotion, floods, epidemics, health and quarantine restrictions, wars, rebellion, transport problems, power supply faults, air-condition faults as well as reasons beyond the control of the company related to supply of manpower, supply of material, manufacturing facilities. In case repairs, setting works, spare parts or maintenance works are required due to the foregoing reasons or due to acts of negligence, default or mistakes by customer or agents thereof.

Variations in costs:

Company is entitled to change costs of maintenance works included in agreement upon completion of contract term or for any following year. Such change shall be notified to customer within (30) days of occurrence thereof. Upon receipt of such notice the customer may

- بالنسبة للإصلاح الذي يحتاج إلى قطع غيار غير متوفرة داخل المملكة فإن الشركة تضمن الإصلاح خلال 8 أسابيع عمل أو أقل من استلام طلب صيانة سواء بالفاكس أو بالإيميل من الشخص المخول بطلب الصيانة من المستشفى .

تضمن الشركة توفير مهندسين وفنيين مدربين بشكل ممتاز على الأجهزة موضوع العقد والمباعدة من الشركة إلى الطرف الأول .

8- نظام تلقي بلاغات الأعطال :-

يقوم الطرف الأول بإبلاغ الطرف الثاني بالأعطال وذلك بواسطة رسالة فاكس أو إيميل وإذا تعذر ذلك يكون البلاغ تليفونيا لقسم الصيانة بالشركة على أن يوضح خلال اتصال الممثل الموجود بالجهاز تاريخ العطل وذلك خلال ساعات العمل اليومي من الساعة الثامنة والنصف صباحا وحتى الواحدة والنصف ظهرا ومن الساعة الثانية والنصف عصرا وحتى الساعة السابعة مساءً عدا يومي الجمعة والسبت من كل أسبوع والعطلات الرسمية .

- رقم خدمة العملاء : 920026426 .

9- قطع الغيار:-

يقوم الطرف الثاني بتوفير قطع الغيار التي يحتاجها الطرف الأول لإصلاح الأجهزة موضوع العقد بناء على تقرير صيانة وعرض أسعار يرفقه مهندس الطرف الثاني ويقبله الطرف الأول

10- خدمات الصيانة:-

سيقوم الطرف الثاني بتوفير خدمات الصيانة المحددة في الاتفاقية للمعدات والأجهزة موضوع العقد وسيتم تقديم هذه الخدمة خلال فترة العمل النهارية العادية خلال الأسبوع المكون من خمسة أيام من الأحد إلى الخميس من الساعة الثامنة والنصف صباحا وحتى الواحدة والنصف ظهرا ومن الساعة الثانية والنصف عصرا وحتى الساعة السابعة مساءً وهي أوقات العمل الرسمية.

و تشمل خدمات الصيانة المقدمة من الشركة ما يلي:-

- القيام بفحوصات خدمات الصيانة المجدولة كما هو محدد من قبل الشركة المصنعة حيث يتم إنجازها حسب الجدول الوارد سابقا في الاتفاقية وتشمل هذه الفحوصات أعمال التشحيم والنظافة وفحوصات وظائف المعدات وفحوصات تسريب التيار وأعمال الضبط والمعايرة ويتم الاتفاق على مواعيد الفحص مسبقا حيث يمنح موظفو الشركة كامل الحرية للوصول إلى المعدات موضوع العقد للقيام بأعمال الفحص وستتوقف أعمال الفحص إذا لم يتمكن موظفو الشركة من الوصول إلى المعدات والأجهزة.

استثناءات:-

لا تشمل هذه الاتفاقية تقديم أعمال الصيانة والإصلاح إذا كان العطل ناتج عن أسباب خارجة عن نطاق قدرة وسيطرة الشركة مثل أمور القضاء والقدر والأعمال الناتجة عن إقفال السلطات المدنية والعسكرية والحروب الأهلية والحرائق الاضطرابات , الفيضانات , الأمراض , قيود الحجر الصحي , الحروب , أعمال الشغب ومشاكل النقل وأعطال في الكهرباء أو عدم الانتظام في تزويدها وأعطال أنظمة التكييف أو أسباب خارجة عن نطاق قدرة الشركة بتوفير العمالة والمواد أو تسهيلات التصنيع وإذا تطلب الأمر إجراء تصليحات أو أي أعمال ضبط أو استبدال القطع أو توفير الصيانة وذلك نتيجة أي من الأسباب الأتفة الذكر أو بسبب إهمال أو تقصير العميل أو نتيجة لارتكاب أي أعمال خاطئة من قبل العميل أو وكلائه .

التغير في التكاليف

للشركة الحق في تغيير تكاليف أعمال الصيانة الواردة في هذه الاتفاقية بعد انتهاء مدة العقد من هذه الاتفاقية أو في أي سنة تالية وذلك بإعطاء العميل إشعار خطي قبل (30) يوما من أحداث ذلك التغيير وللعميل الحق في إنهاء الاتفاقية عند استلام ذلك الإشعار عن طريق إخطار الشركة بذلك قبل تاريخ

terminate the services of the company through a notice to that effect before the effective date of such change. However, in the event the customer fails to terminate the service of the company, the new costs shall be valid effective the date specified by the company on the notice to the customer.

In the event the customer makes any adjustments, changes or additions in equipment included in this agreement, which changes the company deems either to make change in services rendered by the company under this agreement, or affect the performance of the equipment, and then company shall be entitled to terminate this agreement or increase the costs/ rates included hereunder. In case the costs sustained by the company under this agreement should increase after the date of effectiveness due to any change in the law, action or due to internal decrees issued or put in operation after the effective date of this agreement, such extra/ additional costs shall be added to the costs specified in agreement.

Any other conditions of this agreement, including the term thereof, which may be affected by such changes, shall be appropriately amended within any law, jurisdiction or otherwise. In case a petition is submitted to a receiver, bankruptcy or violation to any of the terms of the agreement by the customer, which violation is not duly remedied, then the company is entitled, on own discretion, to declare the customer in default, and accordingly take the following measures:

1. To declare that all amounts are due and payable under the agreement.
2. File claim for application of agreement terms or recover damage due to violation of agreement.

Limits of liabilities and obligations

The company shall in no way be liable for any damage to first party or his patients because of the equipment, whether such damages are special, contingent or partial, due to violation of contract, security or tort (including default, debts and obligations), including the following without limitation: loss of equipment use or any other related equipment, damage of equipment related to basic equipment, loss of capital costs, value of alternative products, facilities, services, power as well as costs of temporary stoppage of equipment and any claims by the customer for such damage.

In the event company provides any consultation or assistance to customer in relation to any equipment covered by the agreement, while such consultation and assistance is not included in agreement, or in case such consultation or assistance is related to equipment not referred to in this agreement, then this consultation or assistance shall not involve the company in any responsibility whether in contract, security or in cause of damage of injury (including negligence, debts and obligations).

Covering of information:

The information, proposals and ideas conveyed from customer to company in relation to performance of services under this contract shall not be considered as secrets or confidential unless expressly proved in writing upon approval of the official representative of the company.

General note

This agreement is considered to be the complete and exclusive statement between the parties and therefore supersedes all prior or concurrent proposals or communications regarding the agreement or any other business not included hereunder shall not be binding on parties. Any assignment, change or amendment in conditions contained herein shall not be binding on company unless rendered in writing and signed by the duly authorized representative of

سريان التغيير و اذا اخفق العميل باعطاء اشعار بانتهاء خدمات الشركة فإن تكاليف الصيانة الجديدة تصبح سارية ابتداء من التاريخ المحدد في إشعار الشركة للعميل.

و في حالة استحداث العمل لأي تعديلات أو تغييرات أو اضافات على الأجهزة والمعدات التي تشملها هذه الاتفاقية و التي من شأنها في نظر الشركة ان تحدث تغيير

في الخدمات التي تقوم بها الشركة ضمن هذه الاتفاقية أو ان تؤثر على أداء الأجهزة والمعدات فيكون للشركة الحق في إنهاء هذه الاتفاقية أو تغيير التعليف الواردة فيها. وفي حالة زيادة التعليف التي تتحملها الشركة حيل تنفيذ هذه الاتفاقية بعد تاريخ سريان مفعولها نتيجة إلى تغيير في القانون أو أي إجراء أو تنظيم أو نتيجة للمراسيم أو القوانين الداخلية التي تصدر أو يجري العمل بها بعد تاريخ ابتداء العمل في الاتفاقية فإنه يجري إضافة التعليف المترتبة إلى التعليف المحددة في الاتفاقية

ان اى شروط أخرى في هذه الاتفاقية بما في ذلك مدة العقد و التي يمكن ان تتأثر بمثل هذه التغييرات سيجري تعديلها بأشكال المناسب ضمن اى قانون او تشريع او غيره في حال تقديم الممثل عربضة لحارس قضائي أو في حالة الإفلاس أو في حالة نقض العميل اى شرط من شروط الاتفاقية و لم يتم معالجة هذا النقض، فيا للشركة الحق حسب اختيارها ان تعلن ان العميل مهممل و مقصر و لها الحق في اتخاذ اى من الخطوات التالية:-

- 1- إعلان كافة مبالغ الإئتمانية المستحقة واجبة الدفع.
- 2- إقامة دعوى لتطبيق بنود الاتفاقية أو الاسترداد ما حصل من أضرار نتيجة لنقض الاتفاقية

حدود المسئوليات و الالتزامات

لن تكون الشركة مسؤولة في أي حال من الأحوال عن أي أضرار يتعرض لها الطرف الأول أو مرضاء بسبب الأجهزة سواء كانت تلك الأضرار خاصة أو عارضة أو جزئية ناجمة عن مخالفة العقد أو الكلفة أو أضرار متعمدة أو غير متعمدة (بما في ذلك التقصير أو الدين أو الالتزامات) لها في ذلك على سبيل المثال لا الحصر الخسارة في استعمال الأجهزة والمعدات وأي معدات ذات علاقة و التلف في المعدات التي على علاقة بالمعدات الأساسية و كلفة راس المال و ثمن المنتجات البديلة و المرافق و الخدمات و الطاقة البديلة و التكاليف الناتجة عن التوقف المؤقت للأجهزة أو أي مطالب أخرى من جانب العميل عن هذه الأضرار

وفي حالة تقديم الشركة للمشورة والمساعدة للعميل فيما يتعلق بأي أجهزة وأردت في الاتفاقية دون أن تكون هذه المساعدة أو المشورة مشمولة بهذه الاتفاقية أو إذا كانت هذه المساعدة تتعلق بأجهزة لم يشار إليها في هذه الاتفاقية، فإن تقديم هذه المشورة والمساعدة لن يحمل الشركة أي مسؤولية سواء كان ذلك في العقد أو القالة أو في حالة الضرر والتلف (بما في ذلك الإسهال أو الدين أو الائتمات) أو غير ذلك.

الإفصاح عن المعلومات

إن يتم اعتبار أي من المعلومات أو اقتراحات أو أي أفكار تنقل من العميل إلى الشركة بخصوص تنفيذ الخدمات ضمن الاتفاقية بمثابة أسرار و لن يتم رفعها إلى الشركة بطريقة سرية ما لم يثبت ذلك و بطريقة خفية وبموافقة ممثل الشركة الرسمي.

ملاحظه عامه

تعتبر هذه الاتفاقية بمثابة التنظيم الكامل والوحيد للعلاقة بين الطرفين ويكون لها الغلبة على أي وثيقة أو عرض أو مراسلات أو أي وعد أو شروط سابقة التاريخ والتي

الاتفاقية والذتي تبقى جميعها غير ملزمة لأي من الطرفين وإن أي تنازل أو تغيير أو تعديل لأي من الشروط الواردة ضمن الاتفاقية لا يعتبر ملزم للشركة ما لم يكن خطيا وموقعا نيابة عن الشركة من جانب أحد ممثليها المعتمدين و الموقعين بذلك

the parties.

The validity of this agreement, its implementation as well as all other matters related thereto, including amendment thereon, shall all be subject to the Laws of the Kingdom of Saudi Arabia.

The demands for service or any part thereof under this contract, which are in satisfaction and performance and not as a further security or obligation of the company, its affiliates, the manufacturer of the complete equipment or any part thereof, shall be in writing before or at the time or providing the service. Failure to notify the company to that effect shall be deemed as a waiver for any similar demands.

The total or partial non validity of an article or a paragraph of this agreement shall not affect the validity or correctness of the remaining articles or paragraphs of the agreement. Any deletion, addition or change in this agreement or any supplement thereof shall promptly render the same as null and void.

This agreement is subject to all related laws and regulations prevailing in the Kingdom of Saudi Arabia.

This contract has been issued in 4 pages, Supplements A and B. each party maintains one signed and stamped copy.

إن سريان مفعول الاتفاقية و تنفيذها وكافة المسائل المتعلقة بها وكذلك اي تعديلات تجري عليها ستكون خاضعة لقوانين المملكة العربية السعودية.

إن تنفيذ الخدمات في ظل هذه الاتفاقية تعتبر إرضاء و تلبية لطلبات العميل ولا يجوز اعتبارها كفالة أو التزام آخر من الشركة أو الشركات التابعة لها أو من مصنع الأجهزة يجب أن تكون طلبات الخدمة خطية قبل أو في الوقت الذي يتم طلب الخدمة فيه. وان عدم إبلاغ الشركة بذلك يعتبر تنازلاً عن طلب الخدمة.

إن عدم سريان المفعول الكلي أو الجزئي لأي مادة أو فقرة في هذه الاتفاقية لن يؤثر على سريان و صحة ما تبقى من مثل تلك المواد أو الفقرات الواردة في هذه الاتفاقية أي حذف أو إضافة أو تغيير في هذه الاتفاقية أو أي من ملحقاتها يؤدي إلى إلغاء هذه الاتفاقية فوراً

ينطبق على هذه الاتفاقية القوانين والأنظمة ذات العلاقة المعمول بها في المملكة العربية السعودية

حرر هذا العقد في 4 صفحات والملحق (أ) و(ب). ومن نسختين أصليتين يحتفظ كل طرف بنسخة موقعة و مختومة.

Success if from Allah

والله ولي التوفيق .

First Party:

Signature:

Stamp:

[Handwritten Signature]



الطرف الأول:

التوقيع:

الختم:

Second Party: Al-Amin Medical Instruments Co. Ltd.

الطرف الثاني: شركة الأمين للتجهيزات الطبية

Signature:

Stamp:

[Handwritten Signature]

التوقيع:

الختم:



MEDICAL EQUIPMENT SUBJECTED TO CONTRACT

Supplement A
MADINAH ROAD CENTER

| Sr. | Equipment | Serial No. | Model | Manufacturer | Location | Dept. | Maintenance Schedule | |
|-----|----------------------------|--------------|--------------|--------------|----------|-------------|----------------------|------------|
| | | | | | | | First PPM | Second PPM |
| 1 | Ophthalmic Unit | NA | S.Executive | Frastema | A-004 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 2 | Patient Examination Chair | NA | S.Executive | Frastema | A-004 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 3 | Indirect Head Light | 2322 | Omega 500 | Heine | A-004 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 4 | Retinoscope | NA | NA | Heine | A-004 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 5 | Direct Ophthalmoscope | NA | NA | Heine | A-004 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 6 | Lensmeter | 502963 | LM600P | Nidek | A-004 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 7 | Chart Projector | 163690 | CP-770 | Nidek | A-004 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 8 | Phoropter/Refractor | 27209 | RT-600 | Nidek | A-004 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 9 | Charger/Indirect light | NA | EN-30 | Heine | A-004 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 10 | Heine Charger/ Retinoscope | NA | NT-200 | Heine | A-004 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 11 | Ophthalmoscope | NA | NA | Heine | A-004 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 12 | Slit Lamp | 900-1157500 | BM 900 | Frastema | A-004 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 13 | Slit Lamp | 900-1177170 | BM 900 | Haag Streit | A-005 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 14 | Ophthalmic Unit | NA | S.Executive | Frastema | A-005 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 15 | Applanation Tonometer | 900-42.92891 | Hagg streit | Frastema | A-005 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 16 | Patient Examination Chair | NA | S.Executive | Frastema | A-005 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 17 | Indirect Head Light | 16288 | Heine | Heine | A-005 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 18 | Retinoscope | NA | NA | Heine | A-005 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 19 | Direct Ophthalmoscope | NA | NA | Heine | A-005 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 20 | Chart Projector | 163690 | CP-770 | Nidek | A-005 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 21 | Phoropter/Refractor | 20985 | RT-600 | Nidek | A-005 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 22 | Heine Charger | NA | EN-30 | Heine | A-005 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 23 | Ophthalmic Unit | NA | M. Polyvisus | Frastema | A-006 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 24 | Slit Lamp | 2120084 | CSO | Frastema | A-006 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 25 | Applanation Tonometer | 900-4292894 | Hagg streit | Frastema | A-006 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 26 | Patient Examination Chair | NA | M. Polyvisus | Frastema | A-006 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 27 | Indirect Head Light | 16295 | Omega 100 | Heine | A-006 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 28 | Retinoscope | NA | NA | Heine | A-006 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 29 | Direct Ophthalmoscope | NA | NA | Heine | A-006 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 30 | Chart Projector | 263059 | CP-670 | Nidek | A-006 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 31 | Phoropter/Refractor | 20954 | RT-600 | Nidek | A-006 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 32 | Charger/Indirect light | 3137 | EN-30 | Heine | A-006 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 33 | Ophthalmic Unit | NA | S.Executive | Frastema | A-007 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 34 | Slit Lamp | 900-1157494 | BM 900 | Haag Streit | A-007 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 35 | Applanation Tonometer | 900-4292885 | Hagg streit | Frastema | A-007 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 36 | Patient Examination Chair | NA | S.Executive | Frastema | A-007 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 37 | Indirect Head Light | 15292 | Omega 100 | Heine | A-007 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 38 | Retinoscope | NA | NA | Heine | A-007 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 39 | Direct Ophthalmoscope | NA | NA | Heine | A-007 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 40 | Exophthalmometer | 180556 | NA | Inani | A-007 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 41 | Chart Projector | 263758 | CP-670 | Nidek | A-007 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 42 | Phoropter/Refractor | 26049 | RT-600 | Nidek | A-007 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 43 | Charger/Indirect light | 2316 | EN-30 | Heine | A-007 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 44 | Heine Charger/ Retinoscope | NA | NT-200 | Heine | A-007 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 45 | Ophthalmoscope | NA | S.Executive | Frastema | A-008 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 46 | Slit Lamp | 900-1177168 | BM 900 | Haag Streit | A-008 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 47 | Applanation Tonometer | 900-42147377 | Hagg streit | Frastema | A-008 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 48 | Patient Examination Chair | NA | S.Executive | Frastema | A-008 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 49 | Indirect Head Light | 1120005583 | Omega 500 | Heine | A-008 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |

| | | | | | | | | |
|-----|----------------------------|--------------|-------------|-------------|--------|-------------|-------------|------------|
| 49 | Retinoscope | NA | NA | Heine | A-008 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 50 | Direct Ophthalmoscope | NA | NA | Heine | A-008 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 51 | Chart Projector | 262160 | CP-670 | Nidek | A-008 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 52 | Phoropter/Refractor | 21858 | RT-600 | Nidek | A-008 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 53 | Heine Charger/ Retinoscope | NA | NT-200 | Heine | A-008 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 54 | Ophthalmoscope | NA | S.Executive | Frastema | A-0035 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 55 | Slit Lamp | 900-11533231 | BM 900 | Haag Streit | A-0035 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 56 | Applanation Tonometer | 42463 | Hagg streit | Frastema | A-0035 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 57 | Patient Examination Chair | NA | S.Executive | Frastema | A-0035 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 58 | Charger Indirect light | 3221 | EN-30 | Heine | A-0035 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 59 | Indirect Head Light | 7762 | Omega 100 | Heine | A-0035 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 60 | Retinoscope | NA | NA | Heine | A-0035 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 61 | Direct Ophthalmoscope | NA | CP-670 | Nidek | A-0035 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 62 | Chart Projector | 262160 | CP-670 | Nidek | A-0035 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 63 | Retinoscope | NA | NA | Heine | A-0035 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 64 | Direct Ophthalmoscope | NA | NA | Heine | A-0035 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 65 | Phoropter/Refractor | 1365 | VT5 | Takagi | A-0035 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 66 | Heine Charger/ Retinoscope | NA | NT-200 | Heine | A-0035 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 67 | Ophthalmoscope | 262824 | CP-670 | Nidek | A-0035 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 68 | Portable slit lamp | 1508706145 | SL-15 | Kowa | A-0035 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 69 | Tonopen XL | 2501357 | Medtronic | Solan | A-0035 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 70 | Charger Indirect light | 3221 | EN-30 | Heine | A-0035 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 71 | Ophthalmic Unit | NA | S.Executive | Frastema | A-0036 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 72 | Slit Lamp | 900-1157501 | BM 900 | Haag Streit | A-0036 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 73 | Applanation Tonometer | AT 900 | Hagg streit | Frastema | A-0036 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 74 | Patient Examination Chair | NA | S.Executive | Frastema | A-0036 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 75 | Indirect Head Light | 12834 | Omega | Heine | A-0036 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 76 | Charger Indirect light | 3315 | EN-30 | Heine | A-0036 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 77 | Chart Projector | 262807 | CP-670 | Nidek | A-0036 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 78 | Retinoscope | NA | NA | Heine | A-0036 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 79 | Direct Ophthalmoscope | NA | NA | Heine | A-0036 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 80 | Heine Charger/ Retinoscope | NA | NT-200 | Heine | A-0036 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 81 | Ophthalmoscope | 113507 | Perkins | Kowa | A-0036 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 82 | Hand Held Tonometer | NA | S.Executive | Frastema | A-0037 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 83 | Slit Lamp | 9001154713 | BM 900 | Haag Streit | A-0037 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 84 | Applanation Tonometer | 9004792886 | Hagg streit | Frastema | A-0037 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 85 | Patient Examination Chair | NA | S.Executive | Frastema | A-0037 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 86 | Indirect Head Light | 13014 | Omega | Heine | A-0037 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 87 | Charger Indirect light | 2886 | EN-30 | Heine | A-0037 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 88 | Phoropter/Refractor | 2369 | VT5 | Takagi | A-0037 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 89 | Chart Projector | 263054 | CP-670 | Nidek | A-0037 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 90 | Retinoscope | NA | NA | Heine | A-0037 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 91 | Direct Ophthalmoscope | NA | NA | Heine | A-0037 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 92 | Heine Charger/ Retinoscope | NA | NT-200 | Heine | A-0037 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 93 | Ophthalmoscope | NA | S.Executive | Frastema | A-0038 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 94 | Slit Lamp | 900-41157033 | BM 900 | Haag Streit | A-0038 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 95 | Applanation Tonometer | 9004292893 | Hagg streit | Frastema | A-0038 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 96 | Patient Examination Chair | NA | S.Executive | Frastema | A-0038 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 97 | Indirect Head Light | 20155 | Omega | Heine | A-0038 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 98 | Phoropter/Refractor | 21050 | RT-600 | Nidek | A-0038 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 99 | Chart Projector | 76221 | CP-690U | Nidek | A-0038 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 100 | Retinoscope | NA | NA | Heine | A-0038 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |

| | | | | | | | | | |
|-----|-------------------------------------|-------------|-------------|----------------|----------------|--------|---------------|-------------|-------------|
| 101 | Direct Ophthalmoscope | NA | NA | NA | Heine | A-0038 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 102 | Heine Charger/ Retinoscope | NA | NT-200 | Heine | Heine | A-0038 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 103 | Ophthalmic Unit | NA | 65 M | Frastema | Frastema | A-0017 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 104 | Slit Lamp | 2100139 | CSO | Haag Streit | Haag Streit | A-0017 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 105 | Patient Examination Chair | NA | NA | Frastema | Frastema | A-0017 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 106 | Retinoscope | NA | NA | Heine | Heine | A-0017 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 107 | Direct Ophthalmoscope | NA | NA | Heine | Heine | A-0017 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 108 | Phoropter/Refractor | 23393 | RT-600 | Nidek | Nidek | A-0017 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 109 | Chart Projector | 760254 | CP-690T | Nidek | Nidek | A-0017 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 110 | Retinoscope | NA | NA | Heine | Heine | A-0017 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 111 | Lensmeter | 502042 | LM600P | Nidek | Nidek | A-0017 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 112 | Keratometer | NA | NA | Magnon | Magnon | A-0017 | OPD-Eye -GF | 07.02.2016. | 5.06.2016. |
| 113 | Autorefractometer | 330802 | ARK 501 A | Nidek | Nidek | A-0027 | Autoref room | 07.02.2016. | 5.06.2016. |
| 114 | Autorefractometer | 334414 | ARK 501 A | Nidek | Nidek | A-0027 | Autoref room | 07.02.2016. | 5.06.2016. |
| 115 | Lensmeter | 502936 | LM-600P | Nidek | Nidek | A-0027 | Autoref room | 07.02.2016. | 5.06.2016. |
| 116 | Lensmeter | 13207 | LM-990 | Nidek | Nidek | A-0027 | Autoref room | 07.02.2016. | 5.06.2016. |
| 117 | Hand Held Autorefractometer | 340736 | ARK 30 | Nidek | Nidek | A-0027 | Autoref room | 07.02.2016. | 5.06.2016. |
| 118 | Hand Held Autorefractometer Printer | 141149 | ARK 30 | Nidek | Nidek | A-0027 | Autoref room | 07.02.2016. | 5.06.2016. |
| 119 | Pachymeter | 2495 | DGH-550 | DGH Technology | DGH Technology | A-0023 | Investigation | 07.02.2016. | 5.06.2016. |
| 120 | Goldman machine | 4035 | 325-UD | Takagi | Takagi | A-0025 | Investigation | 07.02.2016. | 5.06.2016. |
| 121 | Computerized Perimetry | 2148 | Octopus 900 | Hagg Streit | Hagg Streit | A-0025 | Investigation | 09.02.2016. | 07.06.2016. |
| 122 | Automatic reclining chair | 66300 | MP3100D | Takagi | Takagi | A-0026 | Investigation | 07.02.2016. | 5.06.2016. |
| 123 | Echoscan | 10754 | US 3300 | Nidek | Nidek | A-0026 | Investigation | 07.02.2016. | 5.06.2016. |
| 124 | Videographic printer | 56785 | NA | Sony | Sony | A-0026 | Investigation | 07.02.2016. | 5.06.2016. |
| 125 | Ophthalmic Unit | NA | S.Executive | Frastema | Frastema | A-0028 | Pedia 1 | 07.02.2016. | 5.06.2016. |
| 126 | Slit Lamp | 900-117714 | BM 900 | Haag Streit | Haag Streit | A-0028 | Pedia 1 | 07.02.2016. | 5.06.2016. |
| 127 | Applanation Tonometer | AT 900 | Hagg Streit | Frastema | Frastema | A-0028 | Pedia 1 | 07.02.2016. | 5.06.2016. |
| 128 | Patient Examination Chair | NA | S.Executive | Frastema | Frastema | A-0028 | Pedia 1 | 07.02.2016. | 5.06.2016. |
| 129 | Indirect Head Light | 20378 | Omega 500 | Heine | Heine | A-0028 | Pedia 1 | 07.02.2016. | 5.06.2016. |
| 130 | Clear Chart | 400056 | SL-2000 | Nidek | Nidek | A-0028 | Pedia 1 | 07.02.2016. | 5.06.2016. |
| 131 | Retinoscope | NA | NA | Heine | Heine | A-0028 | Pedia 1 | 07.02.2016. | 5.06.2016. |
| 132 | Direct Ophthalmoscope | NA | NA | Heine | Heine | A-0028 | Pedia 1 | 07.02.2016. | 5.06.2016. |
| 133 | Portable slit Lamp | 1508603641 | SL-15 | Kowa | Kowa | A-0028 | Pedia 1 | 07.02.2016. | 5.06.2016. |
| 134 | Hand held Autorefractometer | 340829 | ARK-30 | Nidek | Nidek | A-0028 | Pedia 1 | 07.02.2016. | 5.06.2016. |
| 135 | Hand held printer | 141267 | ARK-30 | Nidek | Nidek | A-0028 | Pedia 1 | 07.02.2016. | 5.06.2016. |
| 136 | Lensmeter | 131676 | LM-1000P | Nidek | Nidek | A-0028 | Pedia 1 | 07.02.2016. | 5.06.2016. |
| 137 | Ophthalmic Unit | NA | S.Executive | Frastema | Frastema | A-0029 | Pedia 2 | 07.02.2016. | 5.06.2016. |
| 138 | Slit Lamp | 7861 | BM 900 | Haag Streit | Haag Streit | A-0029 | Pedia 2 | 07.02.2016. | 5.06.2016. |
| 139 | Applanation Tonometer | 90042138665 | Hagg Streit | Frastema | Frastema | A-0029 | Pedia 2 | 07.02.2016. | 5.06.2016. |
| 140 | Patient Examination Chair | NA | S.Executive | Frastema | Frastema | A-0029 | Pedia 2 | 07.02.2016. | 5.06.2016. |
| 141 | Indirect Head Light | 144722 | Omega 500 | Heine | Heine | A-0029 | Pedia 2 | 07.02.2016. | 5.06.2016. |
| 142 | Clear Chart | 20982-1006 | 13750 | Riechart | Riechart | A-0029 | Pedia 2 | 07.02.2016. | 5.06.2016. |
| 143 | Retinoscope | NA | NA | Heine | Heine | A-0029 | Pedia 2 | 07.02.2016. | 5.06.2016. |
| 144 | Direct Ophthalmoscope | NA | NA | Heine | Heine | A-0029 | Pedia 2 | 07.02.2016. | 5.06.2016. |
| 145 | Portable slit Lamp | 510797 | SL-15 | Kowa | Kowa | A-0029 | Pedia 2 | 07.02.2016. | 5.06.2016. |
| 146 | Hand held Autorefractometer | 340728 | ARK-30 | Nidek | Nidek | A-0029 | Pedia 2 | 07.02.2016. | 5.06.2016. |
| 147 | Hand held printer | 141135 | ARK-30 | Nidek | Nidek | A-0029 | Pedia 2 | 07.02.2016. | 5.06.2016. |
| 148 | Lensmeter | 131748 | LM-1000 | Nidek | Nidek | A-0029 | Pedia 2 | 07.02.2016. | 5.06.2016. |
| 149 | Patient Examination Chair | NA | S.Executive | Frastema | Frastema | A-0030 | Orthophic | 07.02.2016. | 5.06.2016. |
| 150 | Transilluminator light | NA | NA | Weich-Allryn | Weich-Allryn | A-0030 | Orthophic | 07.02.2016. | 5.06.2016. |
| 151 | Charger | NA | NT-300 | Heine | Heine | A-0030 | Orthophic | 07.02.2016. | 5.06.2016. |
| 152 | Visual Acuity Computerized | NA | Machintosh | M&S Tech. | M&S Tech. | A-0030 | Orthophic | 07.02.2016. | 5.06.2016. |

| | | | | | | | | |
|-----|-----------------------------|----------------|--------------|-------------|---------|---------------|-------------|-------------|
| 153 | Visual Acuity Monitor | NA | Smart Screen | M&S Tech. | 2nd flr | Ultrasound rm | 08.02.2016. | 06.06.2016. |
| 154 | Lensmeter | 502380 | LM-600P | Nidek | 2nd flr | Ultrasound rm | 08.02.2016. | 06.06.2016. |
| 155 | Ophthalmic Unit | NA | S.Executive | Frastema | 2nd flr | Retina 1 | 08.02.2016. | 06.06.2016. |
| 156 | Slit Lamp Led Powered | 4055 | BQ900 | Haag Streit | 2nd flr | Retina 1 | 08.02.2016. | 06.06.2016. |
| 157 | Applanation Tonometer | AT 900 | Hagg streit | Frastema | 2nd flr | Retina 1 | 08.02.2016. | 06.06.2016. |
| 158 | Patient Examination Chair | NA | S.Executive | Frastema | 2nd flr | Retina 1 | 08.02.2016. | 06.06.2016. |
| 159 | Indirect Head Light | 1120011030 | Omega 500 | Heine | 2nd flr | Retina 1 | 08.02.2016. | 06.06.2016. |
| 160 | Visual Acuity Monitor | NA | Smart Screen | M&S Tech. | 2nd flr | Retina 1 | 08.02.2016. | 06.06.2016. |
| 161 | w/remote control & keyboard | NA | NA | Heine | 2nd flr | Retina 1 | 08.02.2016. | 06.06.2016. |
| 162 | Retinoscope | NA | NA | Heine | 2nd flr | Retina 1 | 08.02.2016. | 06.06.2016. |
| 163 | Direct Ophthalmoscope | NA | NA | Nidek | 2nd flr | Retina 1 | 08.02.2016. | 06.06.2016. |
| 164 | Refractometer w/ printer | 834784 | RT-5100 | Heine | 2nd flr | Retina 1 | 08.02.2016. | 06.06.2016. |
| 165 | Tonopen | NA | Solan | Riechart | 2nd flr | Retina 1 | 08.02.2016. | 06.06.2016. |
| 166 | Lensmeter | 503110 | LM-600P | Nidek | 2nd flr | Retina 1 | 08.02.2016. | 06.06.2016. |
| 167 | Exophthalmometer | K-0161 | NA | Inani | 2nd flr | Retina 1 | 08.02.2016. | 06.06.2016. |
| 168 | Ophthalmic Unit | NA | Polyvisus | Frastema | 2nd flr | Retina 2 | 08.02.2016. | 06.06.2016. |
| 169 | Slit Lamp | 4053 | BQ900 | Haag Streit | 2nd flr | Retina 2 | 08.02.2016. | 06.06.2016. |
| 170 | Applanation Tonometer | AT 900 | Hagg streit | Frastema | 2nd flr | Retina 2 | 08.02.2016. | 06.06.2016. |
| 171 | Patient Examination Chair | NA | S.Executive | Frastema | 2nd flr | Retina 2 | 08.02.2016. | 06.06.2016. |
| 172 | Indirect Head Light | 1121002883 | Omega 500 | Heine | 2nd flr | Retina 2 | 08.02.2016. | 06.06.2016. |
| 173 | Visual Acuity Monitor | NA | Smart Screen | M&S Tech. | 2nd flr | Retina 2 | 08.02.2016. | 06.06.2016. |
| 174 | w/remote control & keyboard | NA | NA | Heine | 2nd flr | Retina 2 | 08.02.2016. | 06.06.2016. |
| 175 | Retinoscope | NA | NA | Heine | 2nd flr | Retina 2 | 08.02.2016. | 06.06.2016. |
| 176 | Direct Ophthalmoscope | NA | NA | Riechart | 2nd flr | Retina 2 | 08.02.2016. | 06.06.2016. |
| 177 | Tonopen | NA | Solan | Nidek | 2nd flr | Retina 2 | 08.02.2016. | 06.06.2016. |
| 178 | Refractometer | 27683 | RT-600 | Frastema | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 179 | Ophthalmic Unit | NA | Polyvisus | Frastema | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 180 | Slit Lamp | 4053 | BQ900 | Haag Streit | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 181 | Applanation Tonometer | AT 900 | Hagg streit | Frastema | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 182 | Patient Examination Chair | NA | Polyvisus | Frastema | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 183 | Indirect Head Light | 1121002883 | Omega 500 | Heine | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 184 | Visual Acuity Monitor | NA | Smart Screen | M&S Tech. | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 185 | w/remote control & keyboard | NA | NA | Heine | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 186 | Retinoscope | NA | NA | Nidek | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 187 | Refractometer | 27905 | RT-600 | Medtronic | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 188 | Tonopen | NA | XL | Frastema | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 189 | Ophthalmic Unit | NA | Polyvisus | Frastema | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 190 | Slit Lamp | 15001 | BM 900 | Haag Streit | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 191 | Applanation Tonometer | 90042138652 | Hagg streit | Frastema | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 192 | Patient Examination Chair | NA | S.Executive | Frastema | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 193 | Indirect Head Light | 1121002887 | Omega 500 | Heine | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 194 | Visual Acuity Monitor | NA | Smart Screen | M&S Tech. | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 195 | w/remote control & keyboard | NA | NA | Heine | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 196 | Retinoscope | NA | NA | Heine | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 197 | Direct Ophthalmoscope | NA | NA | Medtronic | 2nd flr | Retina 3 | 08.02.2016. | 06.06.2016. |
| 198 | Tonopen | NA | XL | Frastema | 2nd flr | Retina 4 | 08.02.2016. | 06.06.2016. |
| 199 | Ophthalmic Unit | NA | Polyvisus | Frastema | 2nd flr | Retina 4 | 08.02.2016. | 06.06.2016. |
| 200 | Slit Lamp | 15638 | BM 900 | Haag Streit | 2nd flr | Retina 4 | 08.02.2016. | 06.06.2016. |
| 201 | Applanation Tonometer | AT 900 CE 1250 | Hagg streit | Frastema | 2nd flr | Retina 4 | 08.02.2016. | 06.06.2016. |
| 202 | Patient Examination Chair | NA | S.Executive | Frastema | 2nd flr | Retina 4 | 08.02.2016. | 06.06.2016. |
| 203 | Indirect Head Light | 1121000636 | Omega 500 | Heine | 2nd flr | Retina 4 | 08.02.2016. | 06.06.2016. |

| | | | | | | | | |
|-----|---|---------------|---------------|----------------|---------|------------------|-------------|-------------|
| 200 | Visual Acuity Monitor w/remote control & keyboard | NA | Smart Screen | M&S Tech. | 2nd flr | Retina 4 | 08.02.2016. | 06.06.2016. |
| 201 | Refinoscope | NA | NA | Heine | 2nd flr | Retina 4 | 08.02.2016. | 06.06.2016. |
| 202 | Direct Ophthalmoscope | NA | NA | Heine | 2nd flr | Retina 4 | 08.02.2016. | 06.06.2016. |
| 203 | Ophthalmic Unit | NA | Polvisus | Frastema | 2nd flr | Retina 5 | 08.02.2016. | 06.06.2016. |
| 204 | Slit Lamp | 4052 | BQ900 | Haag Streit | 2nd flr | Retina 5 | 08.02.2016. | 06.06.2016. |
| 205 | Applanation Tonometer | AT 900 | Hagg streit | Frastema | 2nd flr | Retina 5 | 08.02.2016. | 06.06.2016. |
| 206 | Patient Examination Chair | NA | Polvisus | Frastema | 2nd flr | Retina 5 | 08.02.2016. | 06.06.2016. |
| 207 | Indirect Head Light | 1121002880 | Omega 500 | Heine | 2nd flr | Retina 5 | 08.02.2016. | 06.06.2016. |
| 208 | Visual Acuity Monitor w/remote control & keyboard | NA | Smart Screen | M&S Tech. | 2nd flr | Retina 5 | 08.02.2016. | 06.06.2016. |
| 209 | Refinoscope | NA | NA | Heine | 2nd flr | Retina 5 | 08.02.2016. | 06.06.2016. |
| 210 | Direct Ophthalmoscope | NA | NA | Heine | 2nd flr | Retina 5 | 08.02.2016. | 06.06.2016. |
| 211 | Refractometer | 27930 | RT-600 | Nidek | 2nd flr | Retina 5 | 08.02.2016. | 06.06.2016. |
| 212 | Ophthalmic Unit | NA | Polvisus | Frastema | 2nd flr | Retina 6 | 08.02.2016. | 06.06.2016. |
| 213 | Slit Lamp | 14998 | BM900 | Haag Streit | 2nd flr | Retina 6 | 08.02.2016. | 06.06.2016. |
| 214 | Applanation Tonometer | CE 0124 | Hagg streit | Frastema | 2nd flr | Retina 6 | 08.02.2016. | 06.06.2016. |
| 215 | Patient Examination Chair | NA | Polvisus | Frastema | 2nd flr | Retina 6 | 08.02.2016. | 06.06.2016. |
| 216 | Indirect Head Light | 1121002884 | Omega 500 | Heine | 2nd flr | Retina 6 | 08.02.2016. | 06.06.2016. |
| 217 | Visual Acuity Monitor w/remote control & keyboard | NA | Smart Screen | M&S Tech. | 2nd flr | Retina 6 | 08.02.2016. | 06.06.2016. |
| 218 | Refinoscope | NA | NA | Heine | 2nd flr | Retina 6 | 08.02.2016. | 06.06.2016. |
| 219 | Direct Ophthalmoscope | NA | NA | Heine | 2nd flr | Retina 6 | 08.02.2016. | 06.06.2016. |
| 220 | Recliner | NA | NA | Steris | 2nd flr | Angiography rm | 08.02.2016. | 06.06.2016. |
| 221 | PDT machine | P0184 | Coherent Opal | Hagg Streit | 2nd flr | OCT room | 08.02.2016. | 06.06.2016. |
| 222 | Slit Lamp | 1841 | BM 900 | Hagg Streit | 2nd flr | OCT room | 08.02.2016. | 06.06.2016. |
| 223 | Argon Krypton Laser machine | 70473A | MC-7000 | Nidek | 2nd flr | Laser room | 09.02.2016. | 07.06.2016. |
| 224 | Slit Lamp | NA | SL 1600 | Hagg Streit | 2nd flr | Laser room | 08.02.2016. | 06.06.2016. |
| 225 | Autorefractometer | 335422 | ARK 501A | Nidek | 2nd flr | N-Assessment | 08.02.2016. | 06.06.2016. |
| 226 | Lensmeter | 503892 | LM-600P | Nidek | 2nd flr | N-Assessment | 08.02.2016. | 06.06.2016. |
| 227 | Visual Acuity Monitor w/remote control & keyboard | NA | Smart Screen | M&S Tech. | 2nd flr | N-Assessment | 08.02.2016. | 07.06.2016. |
| 228 | Multi color Laser Machine | 50160 | MC-500 | Nidek | 2nd flr | Laser Dr. Walid | 09.02.2016. | 07.06.2016. |
| 229 | Slit Lamp | 25008 | SL 1800 | Nidek | 2nd flr | Laser Dr. Walid | 08.02.2016. | 06.06.2016. |
| 230 | Slit Lamp | 6111275 | BM 900 | Nidek | 2nd flr | Laser Dr. Walid | 08.02.2016. | 06.06.2016. |
| 231 | Microscope | 120-6377 | OM-8 | Takagi | 2nd flr | IVT room | 08.02.2016. | 06.06.2016. |
| 232 | Microscope | 09-12377 | OM-8 | Takagi | 2nd flr | IVT room | 08.02.2016. | 06.06.2016. |
| 233 | Excimer laser machine | 560046 | EC-5000 | Quest | 3rd flr | Excimer laser rm | | |
| 234 | Intralase | 0508-60090 | 2 | Nidek | 3rd flr | Excimer laser rm | | |
| 235 | Moria Evolution | 5562 | 19380 | Nidek | 3rd flr | Excimer laser rm | 09.02.2016. | 07.06.2016. |
| 236 | Pachymeter | 2011-3140 | DGH-550 | DGH Technology | 3rd flr | Excimer laser rm | 09.02.2016. | 07.06.2016. |
| 237 | UV-X | 1000-0000-000 | Glob Tek | NA | 3rd flr | Excimer laser rm | 09.02.2016. | 07.06.2016. |
| 238 | Patient bed | 17381 | NA | BKRUS | 3rd flr | Excimer laser rm | | |
| 239 | Patient bed | 122796-24 | MK90/YXZ | Dexta | 3rd flr | Excimer laser rm | | |
| 240 | Medi Pump Aspirator | 9980000124 | 1243 | Thomas | 3rd flr | Excimer laser rm | | |
| 241 | Lensmeter | 31310 | LM-350 | Nidek | 3rd flr | Excimer laser rm | 09.02.2016. | 07.06.2016. |
| 242 | Slit Lamp | 21025 | LM-450 | Hagg Streit | 3rd flr | Excimer laser rm | 09.02.2016. | 07.06.2016. |
| 243 | Slit Lamp | 21028 | LM-450 | Hagg Streit | 3rd flr | Excimer laser rm | 09.02.2016. | 07.06.2016. |
| 244 | Hydraulic stool | 082710-1 | 905 | Dexta | 3rd flr | Excimer laser rm | | |
| 245 | Ophthalmic Unit | NA | Polvisus | Frastema | 3rd flr | Optometrist | 09.02.2016. | 07.06.2016. |
| 246 | Slit Lamp | 0211-12060052 | CSO | Haag Streit | 3rd flr | Optometrist | 09.02.2016. | 07.06.2016. |
| 247 | Patient Examination Chair | NA | Polvisus | Frastema | 3rd flr | Optometrist | 09.02.2016. | 07.06.2016. |

| | | | | | | | |
|------------------------------------|------------|--------------|----------------|---------|---------------|-------------|-------------|
| Visual Acuity Monitor | NA | NA | M&S Tech. | 3rd fir | Optometrist | 09.02.2016. | 07.06.2016. |
| 248 w/remote control & keyboard | NA | NA | Heine | 3rd fir | Optometrist | 09.02.2016. | 07.06.2016. |
| 249 Retinoscope | NA | NA | Heine | 3rd fir | Optometrist | 09.02.2016. | 07.06.2016. |
| 250 Direct Ophthalmoscope | NA | NA | Nidek | 3rd fir | Optometrist | 09.02.2016. | 07.06.2016. |
| 251 Refractometer | 27678 | RT-600 | Nidek | 3rd fir | Optometrist | 09.02.2016. | 07.06.2016. |
| 252 Lensmeter | 503106 | LM-600P | Nidek | 3rd fir | Optometrist | 09.02.2016. | 07.06.2016. |
| 253 Radioscope | 10368 | Magnon | Hagg Streit | 3rd fir | Optometrist | 09.02.2016. | 07.06.2016. |
| 254 Keratometer | 21169 | Magnon | Hagg Streit | 3rd fir | Optometrist | 09.02.2016. | 07.06.2016. |
| 255 Ophthalmic Unit | NA | Polyvisus | Frastema | 3rd fir | Refractive 1 | 09.02.2016. | 07.06.2016. |
| 256 Slit Lamp Led Powered | 1954 | BQ900 | Haag Streit | 3rd fir | Refractive 1 | 09.02.2016. | 07.06.2016. |
| 257 Applanation Tonometer | AT 900 | Hagg streit | Frastema | 3rd fir | Refractive 1 | 09.02.2016. | 07.06.2016. |
| 258 Patient Examination Chair | NA | S.Executive | Frastema | 3rd fir | Refractive 1 | 09.02.2016. | 07.06.2016. |
| 259 Indirect Head Light | 1120009560 | Omega 500 | Heine | 3rd fir | Refractive 1 | 09.02.2016. | 07.06.2016. |
| Visual Acuity Monitor | NA | Smart Screen | M&S Tech. | 3rd fir | Refractive 1 | 09.02.2016. | 07.06.2016. |
| 260 w/remote control & keyboard | NA | NA | Heine | 3rd fir | Refractive 1 | 09.02.2016. | 07.06.2016. |
| 261 Retinoscope | NA | NA | Heine | 3rd fir | Refractive 1 | 09.02.2016. | 07.06.2016. |
| 262 Direct Ophthalmoscope | 833297 | RT-5100 | Nidek | 3rd fir | Refractive 1 | 09.02.2016. | 07.06.2016. |
| 263 Refractometer w/ printer | 502958 | LM-600P | Nidek | 3rd fir | Refractive 1 | 09.02.2016. | 07.06.2016. |
| 264 Lensmeter | 20113210 | DGH-550 | DGH Technology | 3rd fir | Refractive 1 | 09.02.2016. | 07.06.2016. |
| 265 Pachymeter | H02594 | K-0161 | Inani | 3rd fir | Refractive 1 | 09.02.2016. | 07.06.2016. |
| 266 Exophthalmometer | NA | LAMBDA 100 | Nidek | 3rd fir | Refractive 1 | 09.02.2016. | 07.06.2016. |
| 267 Portable PotentialAcuity Meter | NA | Polyvisus | Frastema | 3rd fir | Refractive 2 | 09.02.2016. | 07.06.2016. |
| 268 Ophthalmic Unit | 16214 | BM900 | Haag Streit | 3rd fir | Refractive 2 | 09.02.2016. | 07.06.2016. |
| 269 Slit Lamp | CE 1250 | Hagg streit | Frastema | 3rd fir | Refractive 2 | 09.02.2016. | 07.06.2016. |
| 270 Applanation Tonometer | NA | S.Executive | Frastema | 3rd fir | Refractive 2 | 09.02.2016. | 07.06.2016. |
| 271 Patient Examination Chair | 112001058 | Omega 500 | Heine | 3rd fir | Refractive 2 | 09.02.2016. | 07.06.2016. |
| 272 Indirect Head Light | NA | Smart Screen | M&S Tech. | 3rd fir | Refractive 2 | 09.02.2016. | 07.06.2016. |
| Visual Acuity Monitor | NA | NA | Heine | 3rd fir | Refractive 2 | 09.02.2016. | 07.06.2016. |
| 273 w/remote control & keyboard | NA | NA | Heine | 3rd fir | Refractive 2 | 09.02.2016. | 07.06.2016. |
| 274 Retinoscope | 27672 | RT-600 | Nidek | 3rd fir | Refractive 2 | 09.02.2016. | 07.06.2016. |
| 275 Direct Ophthalmoscope | 503108 | LM-600P | Nidek | 3rd fir | Refractive 2 | 09.02.2016. | 07.06.2016. |
| 276 Refractometer | NA | Polyvisus | Frastema | 3rd fir | Refractive 3 | 09.02.2016. | 07.06.2016. |
| 277 Lensmeter | 16213 | BM900 | Haag Streit | 3rd fir | Refractive 3 | 09.02.2016. | 07.06.2016. |
| 278 Ophthalmic Unit | CE 1250 | Hagg streit | Frastema | 3rd fir | Refractive 3 | 09.02.2016. | 07.06.2016. |
| 279 Slit Lamp | NA | S.Executive | Frastema | 3rd fir | Refractive 3 | 09.02.2016. | 07.06.2016. |
| 280 Applanation Tonometer | 1120009585 | Omega 500 | Heine | 3rd fir | Refractive 3 | 09.02.2016. | 07.06.2016. |
| 281 Patient Examination Chair | NA | Smart Screen | M&S Tech. | 3rd fir | Refractive 3 | 09.02.2016. | 07.06.2016. |
| 282 Indirect Head Light | NA | NA | Heine | 3rd fir | Refractive 3 | 09.02.2016. | 07.06.2016. |
| Visual Acuity Monitor | NA | Smart Screen | M&S Tech. | 3rd fir | Refractive 3 | 09.02.2016. | 07.06.2016. |
| 283 w/remote control & keyboard | NA | NA | Heine | 3rd fir | Refractive 3 | 09.02.2016. | 07.06.2016. |
| 284 Retinoscope | 27682 | RT-600 | Nidek | 3rd fir | Refractive 3 | 09.02.2016. | 07.06.2016. |
| 285 Direct Ophthalmoscope | 503117 | LM-600P | DGH Technology | 3rd fir | Refractive 3 | 09.02.2016. | 07.06.2016. |
| 286 Refractometer | 2233 | DGH-550 | Frastema | 3rd fir | Refractive 3 | 09.02.2016. | 07.06.2016. |
| 287 Lensmeter | 1953 | BQ900 | Haag Streit | 3rd fir | Oculoplasty 1 | 09.02.2016. | 07.06.2016. |
| 288 Pachymeter | CE 1250 | Hagg streit | Frastema | 3rd fir | Oculoplasty 1 | 09.02.2016. | 07.06.2016. |
| 289 Ophthalmic Unit | NA | S.Executive | Frastema | 3rd fir | Oculoplasty 1 | 09.02.2016. | 07.06.2016. |
| 290 Slit Lamp | 112000996 | Omega 500 | Heine | 3rd fir | Oculoplasty 1 | 09.02.2016. | 07.06.2016. |
| 291 Applanation Tonometer | NA | Smart Screen | M&S Tech. | 3rd fir | Oculoplasty 1 | 09.02.2016. | 07.06.2016. |
| 292 Patient Examination Chair | NA | NA | Heine | 3rd fir | Oculoplasty 1 | 09.02.2016. | 07.06.2016. |
| 293 Indirect Head Light | NA | NA | Heine | 3rd fir | Oculoplasty 1 | 09.02.2016. | 07.06.2016. |

| | | | | | | | | |
|-----|---|--------------|----------------|-------------|---------|------------------|-------------|-------------|
| 294 | Visual Acuity Monitor w/remote control & keyboard | NA | Smart Screen | M&S Tech. | 3rd fir | Oculoplasty 1 | 09.02.2016. | 07.06.2016. |
| 295 | Refinoscope | NA | NA | Heine | 3rd fir | Oculoplasty 1 | 09.02.2016. | 07.06.2016. |
| 296 | Direct Ophthalmoscope | NA | NA | Heine | 3rd fir | Oculoplasty 1 | 09.02.2016. | 07.06.2016. |
| 297 | Refractometer w/ printer | 833372 | RT-510 | Nidek | 3rd fir | Oculoplasty 1 | 09.02.2016. | 07.06.2016. |
| 298 | Lensmeter | 503110 | LM-600P | Nidek | 3rd fir | Oculoplasty 1 | 09.02.2016. | 07.06.2016. |
| 299 | Exophthalmometer | 1102595 | NA | Inani | 3rd fir | Oculoplasty 1 | 09.02.2016. | 07.06.2016. |
| 300 | Transilluminator light | NA | NA | Heine | 3rd fir | Oculoplasty 1 | 09.02.2016. | 07.06.2016. |
| 301 | Ophthalmic Unit | NA | Polylus | Frastema | 3rd fir | Oculoplasty 2 | 09.02.2016. | 07.06.2016. |
| 302 | Slit Lamp | 16215 | BM900 | Haag Streit | 3rd fir | Oculoplasty 2 | 09.02.2016. | 07.06.2016. |
| 303 | Applanation Tonometer | CE 1250 | Haga streit | Frastema | 3rd fir | Oculoplasty 2 | 09.02.2016. | 07.06.2016. |
| 304 | Patient Examination Chair | NA | S.Executive | Frastema | 3rd fir | Oculoplasty 2 | 09.02.2016. | 07.06.2016. |
| 305 | Indirect Head Light | 1120009562 | Omega 500 | Heine | 3rd fir | Oculoplasty 2 | 09.02.2016. | 07.06.2016. |
| 306 | Visual Acuity Monitor w/remote control & keyboard | NA | Smart Screen | M&S Tech. | 3rd fir | Oculoplasty 2 | 09.02.2016. | 07.06.2016. |
| 307 | Refinoscope | NA | NA | Heine | 3rd fir | Oculoplasty 2 | 09.02.2016. | 07.06.2016. |
| 308 | Direct Ophthalmoscope | NA | NA | Heine | 3rd fir | Oculoplasty 2 | 09.02.2016. | 07.06.2016. |
| 309 | Refractometer | 27683 | RT-500 | Nidek | 3rd fir | Oculoplasty 2 | 09.02.2016. | 07.06.2016. |
| 310 | Lensmeter | 502963 | LM-600P | Nidek | 3rd fir | Oculoplasty 2 | 09.02.2016. | 07.06.2016. |
| 311 | Exophthalmometer | 174090 | NA | Inani | 3rd fir | Oculoplasty 2 | 09.02.2016. | 07.06.2016. |
| 312 | Microscope | 1199538 | OM-10 | Takagi | 3rd fir | Minor room | 09.02.2016. | 07.06.2016. |
| 313 | Surgical floor lamp | 265 | Million lights | Tokyo Atago | 3rd fir | Minor room | 09.02.2016. | 07.06.2016. |
| 314 | Recliner | 73E0000AUB | ESDI-ST | Steris | 3rd fir | Minor room | 09.02.2016. | 07.06.2016. |
| 315 | Headlight | 1/1319 | HK-7000 | Heini | 3rd fir | Minor room | 09.02.2016. | 07.06.2016. |
| 316 | Hyfrecator | 13EGJ364 | 2000 | Commed | 3rd fir | Minor room | 09.02.2016. | 07.06.2016. |
| 317 | Electrolysis | 9222 | B2122 | John Weiss | 3rd fir | Minor room | 09.02.2016. | 07.06.2016. |
| 318 | Dermo Laser Machine w/ dioxide cylinder goggles & carboxy cream | 13CTR302-875 | Carboxy | Won-Cosjet | 3rd fir | Dermo room | 09.02.2016. | 07.06.2016. |
| 319 | Photo Slit Lamp w/ Biomicroscope /Canon Camera flash unit/Slit lamp illuminator | 388 | BX 900 | Hagg Streit | 3rd fir | Investigation rm | 09.02.2016. | 07.06.2016. |
| 320 | Scanning Endothelial System w/ laptop/HP printer | 4571 | SEM Technology | Nidek | 3rd fir | Investigation rm | 09.02.2016. | 07.06.2016. |
| 321 | Cornel Endothelial machine w/ laptop/HP printer | 140406 | CEM-530 | Nidek | 3rd fir | Investigation rm | 09.02.2016. | 07.06.2016. |
| 322 | Cornel Topography | 30218 | OPD Scan II | Nidek | 3rd fir | Investigation rm | 09.02.2016. | 07.06.2016. |
| 323 | Autorefractometer | 334925 | ARK 510A | Nidek | 3rd fir | N.Assessment rm | 09.02.2016. | 07.06.2016. |
| 324 | Portable suction machine | 173449/ASP | NA | MGM | 3rd fir | N.Assessment rm | 09.02.2016. | 07.06.2016. |

| Equipment | Examination Unit | Check list |
|--------------------|------------------|------------|
| Date | | |
| Location | | |
| Service Engineer | | |
| Customer | | |
| Main body | SN: | |
| SL- | SN: | |
| RT- | SN: | |
| Chart Projector CP | SN: | |
| Indirect | SN: | |

- ☐ First PPM
☐ Second PPM

| Special Tools Needed | Calibration Date |
|------------------------|------------------|
| Special Cleaning Paper | |
| Methanol | |
| | |
| | |
| | |
| | |
| | |
| | |

- ☐ Checking Indirect Bulb
☐ Attaching a pad of the chin rest
☐ Checking the fuses
☐ Checking the illumination
☐ Cleaning outer parts
☐ Cleaning the optics
☐ Checking the chair movement
☐ Checking the head light
☐ Checking the chart projector illumination
☐ Checking the phoropter
☐ Checking the table movement

Remark

| |
|--|
| |
| |

| | |
|---------------------------|--|
| Customer Name: | |
| Customer/ Dr. Signature : | |

| | |
|----------------------|-------------------------------|
| Equipment | OCTOPUS 900 Check list |
| Date | |
| Location | |
| Service Engineer | |
| Customer | |
| OCTOPUS900 Main Body | |
| Accessory1 | |
| | |

- ☐ First PPM
☐ Second PPM

| | |
|------------------------|------------------|
| Special Tools Needed | Calibration Date |
| Special Cleaning Paper | |
| Methanol | |
| | |
| | |
| | |
| | |
| | |

- ☐ Check Focusing of Target reticle (circle)
☐ Check Characters on the monitor
☐ Wipe surface with alcohol
☐ Clean Dome
☐ Check measurement
☐ Check Proper measurement of Left & Right
☐ Check Chin Rest vertical Smooth movement
☐ Check communication & PC Capacity
☐ Cleaning of optics.

Remark

| |
|--|
| |
| |

| | |
|----------------|--|
| Customer Name: | |
|----------------|--|

| | |
|---------------------------|--|
| Customer/ Dr. Signature : | |
|---------------------------|--|

PM Check List NIDEK CEM Specular Microscope

Equipment:

Customer / Location:

Model :

S/N:

| | Pass | Fail |
|--|--------------------------|--------------------------|
| 1. Focusing of Target reticle (circle) | <input type="checkbox"/> | <input type="checkbox"/> |
| 2. Characters on the T.V monitor | <input type="checkbox"/> | <input type="checkbox"/> |
| 3. Smooth horizontal movement of main body | <input type="checkbox"/> | <input type="checkbox"/> |
| 4. Clean the lenses | <input type="checkbox"/> | <input type="checkbox"/> |
| 5. Test real eye | <input type="checkbox"/> | <input type="checkbox"/> |
| 6. Proper measurement of spectrometer | <input type="checkbox"/> | <input type="checkbox"/> |
| 7. Chin Rest vertical Smooth movement | <input type="checkbox"/> | <input type="checkbox"/> |
| 8. Cleaning of optics. | <input type="checkbox"/> | <input type="checkbox"/> |

Unit meets all specifications and suitable for its intended use. () Yes () No

Test performed by _____
Date: _____

Signature _____

Notes:

Equipment **MC500 VIXI Check list**

| | |
|------------------|--|
| Date | |
| Location | |
| Service Engineer | |
| Customer | |
| Slit lamp | |
| Accessory1 | |
| Accessory2 | |

- ☐ First PPM
☐ Second PPM

| | |
|-----------------------|------------------|
| Special Tools Needed | Calibration Date |
| Power meter | |
| Methanol | |
| Special Clining Paper | |
| | |
| | |
| | |
| | |

Power Data

| Color | Power indicator | Spot size | |
|------------------|-----------------|-----------|-------|
| GREEN,YELLOW,RED | | 50µm | 500µm |
| | 50 mw | | |
| | 100 mw | | |
| | 200 mw | | |
| | 500 mw | | |
| | 900 mw | | |
| | 1000 mw | | |

- ☐ Checking the fiber alignment
☐ Checking the indicated spot size
☐ Checking the protective filter
☐ Checking the laser output
☐ Cleanning the optics
☐ Checking the output and Current of LH RED -----A Green -----A Yellow -----A
☐ Checking the power source

Remark

| |
|--|
| |
| |

Customer Name:

Customer/ Dr. Signature :

MC-7000 PPM CHECK LIST

| Customer | Magrabi Hospital |
|-------------------|------------------|
| MC-7000 main body | |
| Slit lamp | S/N |
| Accessory1 | S/N |
| Accessory2 | S/N |

| Date | |
|------------------|--|
| Location | |
| Service Engineer | |

spot size : 500µm

Power Data

| Color | Power indicator | Power Meter | | |
|--------------------|-----------------|-------------|-----|-----|
| | | RED | YEL | GYN |
| GREEN, YELLOW, RED | | | | |
| | 50 mw | | | |
| | 100 mw | | | |
| | 200 mw | | | |
| | 500 mw | | | |
| | 900 mw | | | |
| | 1000 mw | | | |

- ☐ Checking the fiber alignment
- ☐ Checking the indicated spot size
- ☐ Checking the protective filter
- ☐ Checking the laser output
- ☐ Cleaning the optics
- ☐ Checking the output and volatage of LH
- ☐ Checking the power source
- ☐ Checking the water level
- ☐ Checking the footswitch

Remarks

US-3300 PPM CHECK LIST

| Customer | |
|-------------------|------|
| US-3300 main body | S/N: |
| A PROBE | S/N |
| B PROBE | S/N |
| Accessory | S/N |

| | |
|------------------|--|
| Date | |
| Location | |
| Service Engineer | |

- ☐ Measure the A Probe reading by phantom crystal /
- ☐ Gain control knob Rotation
- ☐ Check B probe & cable
- ☐ Check freezing in A&B-mode
- ☐ Check footswitch
- ☐ Checking noise in B mode
- ☐ Checking the power source

Remarks

DGH Pachette Checklist (calibration)

| | |
|-----------|-----|
| Customer | |
| Main body | SN: |

| |
|-------------------|
| Date: |
| Location: |
| Service Engineer: |

- ☐ Probe connected to front panel
- ☐ Probe is connected to front panel
- ☐ Power switch is in the Off ("0") position
- ☐ While pressing in the Up arrow switch turn on Unit
- ☐ Verify double beep occurs and display indicates as per Calibration Box
- ☐ Disconnect probe and connect calibration box to unit
- ☐ Hold down push button on calibration box
- ☐ Press down on measure side of foot pedal
- ☐ Observe Actual measurement on display (800 +/- 3 nm)
- ☐ Release push button
- ☐ Disconnect calibration box and connect probe. Press Clear

Remarks

| |
|--|
| |
| |
| |
| |

OPD Scan Checklist

| | |
|-----------|-----|
| Customer | |
| Main body | SN: |

| |
|-------------------|
| Date: |
| Location: |
| Service Engineer: |

- ☐ Wipe the measuring window lightly with alcohol mixed solution
- ☐ AR calibration
- ☐ perform the calibrations in 2D increments from -20D to +22D
- ☐ Cleaning outer parts
- ☐ Checking Joystick Movement
- ☐ Perform the AR measurement
- ☐ check that the measurement values are within specifications.
- ☐ Check the mapping function
- ☐ CT Calibration
- ☐ Measure the model eyes from R5.00 to R10.00

| Corneal curvature radius (R) | Tolerance (mm) | Difference between R1 and R2 (mm) |
|------------------------------|----------------|-----------------------------------|
| R 5.00 mm | ±0.02 | 0.04 |
| R 6.50 mm | ±0.02 | 0.04 |
| R 8.00 mm | ±0.02 | 0.04 |
| R 9.40 mm | ±0.03 | 0.04 |
| R 10.00 mm | ±0.04 | 0.04 |

Remarks

.....

.....